

*Medilime*<sup>®</sup>



**Product Catalogue**

## Our vision

MediLime focuses on innovating and improving enteral feeding systems, to make the life of users more enjoyable. MediLime not only co-operates with specialists, such as doctors and nurses, but also with end-users. This enables MediLime to continuously optimize and expand its range.



# Nutrition is essential

Not only to keep our body healthy, but also to enable us to grow and heal. However, there are various medical reasons which can make consuming food difficult. For this group of people, MediLime has developed new products.

**“ I’m very enthusiastic! The button and sealing caps are easier and faster to insert and more pleasant to replace, I see much less irritation.”**

Heleen Voogt is a passionate nursing consultant at the Pediatric Stomach-Liver-Intestine Department at the Academic Medical Center (AMC) in Amsterdam.



## The products

MediLime supplies all types of products for enteral feeding, such as Nasogastric Feeding Tubes, Tummy Buttons, Gastric Feeding Tubes and Feeding Syringes. In addition, our range also has accessories such as caps, clips and connectors.

The MediLime Tummy Button and all associated accessories meet the new global standards for Feeding Catheters. These new products – in contrast to other products – contain only tissue-friendly materials and are free of plasticizers.

## Accessories

MediLime supplies a large number of (ENFit) accessories. The range includes Feeding Syringes, bottle adapters, connectors, sealing caps and clips.



## Extension Tubes

MediLime Extension Tubes are used to administer enteral nutrition and/or medication through the button.

The MediLime Extension Tubes can be used for all buttons and are available with ENFit, ENFit/Y and bolus connectors.

The MediLime Extension Tubes are made of the highest quality polyurethane, are latex free and contain no plasticizers.

The Extension Tubes are available in 30 and 60 cm lengths (12 and 24 inches).

## Connections

The same applies to the Tummy Button, the feedback from end-users and specialists is leading in innovating the Extension Tubes. In order to meet the demand in different settings there is a variety of possible connections:




A photograph of a family of three in a hospital bed. A woman with curly hair is leaning over a young child who is lying in the bed, smiling. A man is sitting next to the child, also smiling. The scene is warm and supportive.

## Table of contents

- 6 MediLime Tummy Button
- 9 MediLime Extension Tubes
- 10 MediLime Stoma Measuring Catheter
- 11 MediLime Gastric Feeding Tubes
- 12 MediLime Jejunal Catheters
- 13 MediLime ENFit Feeding Syringe
- 14 Accessories
  - 15 Bottle adapters
  - 15 ENFit 3-way stopcock
  - 15 Sealing caps
  - 16 Connectors
  - 16 Clips
  - 17 Stoma plugs
  - 17 Drainage bags
  - 17 PH strips
- 18 MediLime EnClean brush

“  
Our ultimate goal  
is to provide  
the best product  
for everyone.”

“  
Making  
life easier.”



## The MediLime Tummy Button

A button can be a comfortable solution, because of its small size and the fact that you don't have to carry a tube. We are committed to listen to the end-users and specialists and we show this by innovating our products year in, year out.

### Some key features

**The MediLime Tummy Button has one sole purpose, and that is to make living with a feeding system more pleasant and efficient.**

- The MediLime Tummy Button consists of pure silicone, is latex- and DEHP free, resulting in a healthy skin around the gastro stoma
- An apple shaped balloon improves an ideal shape for closing the stoma channel
- Best compatibility with all MediLime Extension Tubes
- Various balloon diameters (per size 12FR up to 24FR, from 3 up to 8ml volume) and is available in all lengths (0.8, 1.0, 1.2, 1.5, 1.7, 2.0, 2.3, 2.5, 2.7, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5 and 6.0)
- The Tummy Button comes with an input guide. This can be a useful tool when placing the button if there is scar tissue at the fistula opening or with a larger size of the Tummy Button

### Package contents

The Tummy Button is packed in a box along with the following accessories for easy placement:

- Input guide
- Syringe
- Lubricating gel 5 grams
- Gauze





Article	Item description
<b>12 french</b>	<b>3ml balloon</b>
0320-12-0.8-02	MediLime Tummy Button 12fr 0,8cm
0320-12-1.0-02	MediLime Tummy Button 12fr 1,0cm
0320-12-1.2-02	MediLime Tummy Button 12fr 1,2cm
0320-12-1.5-02	MediLime Tummy Button 12fr 1,5cm
0320-12-1.7-02	MediLime Tummy Button 12fr 1,7cm
0320-12-2.0-02	MediLime Tummy Button 12fr 2,0cm
0320-12-2.3-02	MediLime Tummy Button 12fr 2,3cm
0320-12-2.5-02	MediLime Tummy Button 12fr 2,5cm
0320-12-2.7-02	MediLime Tummy Button 12fr 2,7cm
0320-12-3.0-02	MediLime Tummy Button 12fr 3,0cm
0320-12-3.5-02	MediLime Tummy Button 12fr 3,5cm
0320-12-4.0-02	MediLime Tummy Button 12fr 4,0cm

Article	Item description
<b>14 french</b>	<b>4ml balloon</b>
0320-14-0.8-02	MediLime Tummy Button 14fr 0,8cm
0320-14-1.0-02	MediLime Tummy Button 14fr 1,0cm
0320-14-1.2-02	MediLime Tummy Button 14fr 1,2cm
0320-14-1.5-02	MediLime Tummy Button 14fr 1,5cm
0320-14-1.7-02	MediLime Tummy Button 14fr 1,7cm
0320-14-2.0-02	MediLime Tummy Button 14fr 2,0cm
0320-14-2.3-02	MediLime Tummy Button 14fr 2,3cm
0320-14-2.5-02	MediLime Tummy Button 14fr 2,5cm
0320-14-2.7-02	MediLime Tummy Button 14fr 2,7cm
0320-14-3.0-02	MediLime Tummy Button 14fr 3,0cm
0320-14-3.5-02	MediLime Tummy Button 14fr 3,5cm
0320-14-4.0-02	MediLime Tummy Button 14fr 4,0cm
0320-14-4.5-02	MediLime Tummy Button 14fr 4,5cm
0320-14-5.0-02	MediLime Tummy Button 14fr 5,0cm
0320-14-5.5-02	MediLime Tummy Button 14fr 5,5cm
0320-14-6.0-02	MediLime Tummy Button 14fr 6,0cm

Article	Item description
<b>16 french</b>	<b>5ml balloon</b>
0320-16-0.8-02	MediLime Tummy Button 16fr 0,8cm
0320-16-1.0-02	MediLime Tummy Button 16fr 1,0cm
0320-16-1.2-02	MediLime Tummy Button 16fr 1,2cm
0320-16-1.5-02	MediLime Tummy Button 16fr 1,5cm
0320-16-1.7-02	MediLime Tummy Button 16fr 1,7cm
0320-16-2.0-02	MediLime Tummy Button 16fr 2,0cm
0320-16-2.3-02	MediLime Tummy Button 16fr 2,3cm
0320-16-2.5-02	MediLime Tummy Button 16fr 2,5cm
0320-16-2.7-02	MediLime Tummy Button 16fr 2,7cm
0320-16-3.0-02	MediLime Tummy Button 16fr 3,0cm
0320-16-3.5-02	MediLime Tummy Button 16fr 3,5cm
0320-16-4.0-02	MediLime Tummy Button 16fr 4,0cm
0320-16-4.5-02	MediLime Tummy Button 16fr 4,5cm
0320-16-5.0-02	MediLime Tummy Button 16fr 5,0cm

Article	Item description
<b>18 french</b>	<b>5ml balloon</b>
0320-18-1.0-02	MediLime Tummy Button 18fr 1,0cm
0320-18-1.2-02	MediLime Tummy Button 18fr 1,2cm
0320-18-1.5-02	MediLime Tummy Button 18fr 1,5cm
0320-18-1.7-02	MediLime Tummy Button 18fr 1,7cm
0320-18-2.0-02	MediLime Tummy Button 18fr 2,0cm
0320-18-2.3-02	MediLime Tummy Button 18fr 2,3cm
0320-18-2.5-02	MediLime Tummy Button 18fr 2,5cm
0320-18-2.7-02	MediLime Tummy Button 18fr 2,7cm
0320-18-3.0-02	MediLime Tummy Button 18fr 3,0cm
0320-18-3.5-02	MediLime Tummy Button 18fr 3,5cm
0320-18-4.0-02	MediLime Tummy Button 18fr 4,0cm
0320-18-4.5-02	MediLime Tummy Button 18fr 4,5cm
0320-18-5.0-02	MediLime Tummy Button 18fr 5,0cm

Article	Item description
<b>20 french</b>	<b>8ml balloon</b>
0320-20-0.8-02	MediLime Tummy Button 20fr 0,8cm
0320-20-1.0-02	MediLime Tummy Button 20fr 1,0cm
0320-20-1.2-02	MediLime Tummy Button 20fr 1,2cm
0320-20-1.5-02	MediLime Tummy Button 20fr 1,5cm
0320-20-1.7-02	MediLime Tummy Button 20fr 1,7cm
0320-20-2.0-02	MediLime Tummy Button 20fr 2,0cm
0320-20-2.3-02	MediLime Tummy Button 20fr 2,3cm
0320-20-2.5-02	MediLime Tummy Button 20fr 2,5cm
0320-20-2.7-02	MediLime Tummy Button 20fr 2,7cm
0320-20-3.0-02	MediLime Tummy Button 20fr 3,0cm
0320-20-3.5-02	MediLime Tummy Button 20fr 3,5cm
0320-20-4.0-02	MediLime Tummy Button 20fr 4,0cm
0320-20-4.5-02	MediLime Tummy Button 20fr 4,5cm
0320-20-5.0-02	MediLime Tummy Button 20fr 5,0cm

Article	Item description
<b>24 french</b>	<b>8ml balloon</b>
0320-24-1.5-02	MediLime Tummy Button 24fr 1,5cm
0320-24-1.7-02	MediLime Tummy Button 24fr 1,7cm
0320-24-2.0-02	MediLime Tummy Button 24fr 2,0cm
0320-24-2.3-02	MediLime Tummy Button 24fr 2,3cm
0320-24-2.5-02	MediLime Tummy Button 24fr 2,5cm
0320-24-2.7-02	MediLime Tummy Button 24fr 2,7cm
0320-24-3.0-02	MediLime Tummy Button 24fr 3,0cm
0320-24-3.5-02	MediLime Tummy Button 24fr 3,5cm
0320-24-4.0-02	MediLime Tummy Button 24fr 4,0cm
0320-24-4.5-02	MediLime Tummy Button 24fr 4,5cm
0320-24-5.0-02	MediLime Tummy Button 24fr 5,0cm

## Product features

Material	Silicone
DEHP free	YES
Latex free	YES
Disposable	YES
Sterile	YES
Sterilization method	EO
Visibility	RX
Packaging	Box of blister pack
Type barcode	GS1
Unit	Per piece

# The MediLime Extension Tubes

Our Extension Tubes are used to administer enteral nutrition and/or medication through the button. The MediLime Extension Tubes can be used for all buttons and are available with ENFit, ENFit/Y and bolus connectors.

The MediLime Extension Tubes are made of the highest quality polyurethane, are latex free and contain no plasticizers. The Extension Tubes are available in the lengths 30 cm and 60 cm (12 inches and 24 inches).



- High flow button connection (ideal for usage with a feeding pump)
- Improved material to reduce the impact of medication
- MediLime is the only supplier of Extension Tubes with the purple connections
- Perfectly applicable on all buttons
- More grip for easy placement

Article	Description	Quantity
0328-12	MediLime Extension Tube	1. ENFit
0328-24	MediLime Extension Tube	1. ENFit
0329-12	MediLime Extension Tube	2. Y-ENFit
0329-24	MediLime Extension Tube	2. Y-ENFit
0324-12	MediLime Extension Tube	3. Large ENFit
0324-24	MediLime Extension Tube	3. Large ENFit
0323-12	MediLime Extension Tube	4. Bolus
0323-24	MediLime Extension Tube	4. Bolus

## Connections

Just like the Tummy Button, the feedback from end-users and specialists is leading in innovating the Extension Tubes. In order to meet the demand in different settings there is a variety of possible connections:



## Swivel connector

MediLime Extension Tubes are equipped with a swivel connector. The connector rotates separately from the tube, allowing it to be connected and disconnected more easily.

## Product features

Material	Polyurethane
PHT and BPA free	YES
Latex free	YES
Disposable	YES
Sterile	NO
Packaging	Box with 10 separate packages
Type barcode	GS1
Unit	Per 10 pieces

# MediLime Stoma Measuring Catheter

The MediLime Stoma Measuring Catheter is used to obtain the correct size of the MediLime Tummy Button. The stoma canal is measured with this catheter.

The use of the Stoma Measuring Catheter will result in a more precise Tummy Button placement. We strongly recommend using one for each (new) placement of a button in order to optimize the perfect fit.

Article	Description	Quantity
0325-00	MediLime Stoma Measuring Catheter	incl. syringe 2 pieces



“ Our product developers develop the new products together with specialists and end-users.”

## Product features

Material	Silicone
PHT and BPA free	YES
Latex free	YES
Disposable	YES
Sterile	YES
Packaging	Packed per unit
Type barcode	GS1
Unit	Per 2 pieces

# MediLime Gastric Feeding Tubes

The MediLime Gastric Feeding Tube with ENFit connection is made for administering enteral nutrition and medication. The Gastric Feeding Tube is made of pure silicone and is therefore skin-friendly and pliable. The size of the tube and the balloon capacity are clearly stated on the shaft.



Article	Description	Quantity
0503-12	MediLime Gastric Feeding Tube	1-way with ENFit adapter 12FR 2 pieces
0503-14	MediLime Gastric Feeding Tube	1-way with ENFit adapter 14FR 2 pieces
0503-16	MediLime Gastric Feeding Tube	1-way with ENFit adapter 16FR 2 pieces
0503-18	MediLime Gastric Feeding Tube	1-way with ENFit adapter 18FR 2 pieces
0503-20	MediLime Gastric Feeding Tube	1-way with ENFit adapter 20FR 2 pieces
0503-22	MediLime Gastric Feeding Tube	1-way with ENFit adapter 22FR 2 pieces
0503-24	MediLime Gastric Feeding Tube	1-way with ENFit adapter 24FR 2 pieces



## Product features

Material	Silicone
PHT and BPA free	YES
Latex free	YES
Disposable	YES
Sterile	YES
Packaging	Packed per unit
Type barcode	GS1
Unit	Per 2 pieces

# MediLime Jejunal Feeding Tubes

The MediLime Jejunal Feeding Tube is a 100% silicone device for delivery of enteral nutrition into the distal duodenum or jejunum. The device has a double lumens tube, an air-vented skin disk and a fixation balloon.



Article	Description	Quantity
103080031280A	MediLime Jejunal Feeding Tube 1-way 12FR 5ml	2 pieces
103080031480A	MediLime Jejunal Feeding Tube 1-way 14FR 5ml	2 pieces
103080031680A	MediLime Jejunal Feeding Tube 1-way 16FR 10ml	2 pieces

## Product features

Material	Silicone
PHT and BPA free	YES
Latex free	YES
Disposable	YES
Sterile	YES
Packaging	Packed per unit
Type barcode	GS1
Unit	Per 2 pieces

# MediLime ENFit Feeding Syringe

MediLime offers a complete range of Feeding Syringes with an ENFit connection. These Feeding Syringes are easy to connect to the other ENFit products. It is advised to use MediLime bottle adapters for proper dosing.

Article	Description	Quantity
0909-010	MediLime ENFit Feeding Syringe 1 ml, low dose, with wings	100 pieces
0909-030	MediLime ENFit Feeding Syringe 3ml, low dose, with wings	100 pieces
0909-050	MediLime ENFit Feeding Syringe 5ml, low dose	100 pieces
0909-10	MediLime ENFit Feeding Syringe 10ml	100 pieces
0909-20	MediLime ENFit Feeding Syringe 20ml	50 pieces
0909-60	MediLime ENFit Feeding Syringe 60ml	30 pieces



## Product features

Fitting	ENFit
DEHP free	YES
Latex free	YES
Disposable	YES
Sterile	YES
Packaging	Blister packs
Unit	Differs per size



## Accessories

MediLime supplies a large number of (ENFit) accessories for the administration of enteral nutrition and medication through the Feeding Tube. The range includes bottle adapters, connectors, sealing caps and clips.

MediLime is continuously working on optimizing its range and strives to bring the best product for everyone.

## Bottle adapters

To extract medication from bottles, a bottle adapter is used. This is a screw cap that is placed in the medication bottle which makes it possible to suck the medication with the help of an Feeding Syringe.

The MediLime bottle adapters are available in various diameters for the various sizes of medication bottles.

Article	Description	Quantity
 BAISO1	MediLime ENFit bottle adapter 1 ø 12mm - 15mm	100 pieces
 BAISO2	MediLime ENFit bottle adapter 2 ø 14mm - 16.5mm	100 pieces
 BAISO3	MediLime ENFit bottle adapter 3 ø 17mm - 20.5mm	100 pieces
 BAISO4	MediLime ENFit bottle adapter 4 ø 24mm - 26mm	100 pieces
 SBAISO	MediLime ENFit Stepped bottle adapter, ø 15mm - 27mm	50 pieces

## ENFit 3-way stopcock

The MediLime ENFit 3-way stopcock allows you to connect up to three ENFit products at once. Made from tritan, which makes it more resistant to medication.

Article	Description	Quantity
0800-20	MediLime ENFit stopcock	20 pieces



## Sealing caps

MediLime supplies various sealing caps for both Feeding Syringes and Extension Tubes.

Article	Description	Quantity
0910-01	MediLime sealing cap Sterile, packed per piece	100 pieces
0910-10	MediLime sealing cap Sterile, packed per 10	100 pieces
C00579	MediLime ENFit female cap	30 pieces













## Connectors

MediLime has a wide range of connectors in its product range. With a connector, products with different fittings (ENFit, luer lock, bolus) can be linked together.

Article	Description	Quantity
 10999382	MediLime ENFit female to luerlock female connector	10 pieces
 10999383	MediLime ENFit male to luerlock female connector	10 pieces
 C00539	MediLime ENFit Male to bolus connector	5 pieces
 F00106	MediLime ENFit female to bolus connector	1 piece
 F00106-30	MediLime ENFit female to bolus connector	30 pieces

## Clips

MediLime has various types of clips in its product range. There are clips that are clamped around the tube of the Feeding Tube with or without side notches and separate clips that are clamped around the tube.

Article	Description	Quantity
 13020	MediLime clip for tubes up to 20 FR	10 pieces



## Stoma plugs

MediLime silicone stoma plug sets are designed for use in emergencies to prevent stomas from closing when a gastrostomy catheter or button has fallen out.

Article	Description	Quantity
ENP-S	MediLime silicone stoma plug set 4cm 10,12,14 & 16FR	1
ENP-L	MediLime silicone stoma plug set 7cm 10,12,14 & 16FR	1
ENP-XL	MediLime silicone stoma plug set 7cm / 18 & 20 FR	1

## Drainage bags



MediLime drainage bags with ENFit connectors are suitable for siphoning the stomach. The bags are made of soft and pliable material and are available in various sizes.

Article	Description	Quantity
DBISO-50	MediLime drainage bags 50ml ENFit	20
DBISO-250	MediLime drainage bags 250ml ENFit	20
DBISO-500	MediLime drainage bags 500ml ENFit	20
DBISO-1000	MediLime drainage bags 1000ml ENFit	20

## PH strips



The MediLime PH indicator strips are used to measure the acidity of released fluids when placing a button or Nasogastric Feeding Tube.



Article	Description	Quantity
ENTPH	MediLime PH indicator strips (pH 2-9)	100
PHX0-6	MediLime PH indicator strips (pH 0-6)	100

# MediLime cleaning brushes

The cleaning brushes help you clean out ENFit connectors when there is residue of food or medication. MediLime has a reusable and a disposable brush in its product range. The reusable brush (MediLime EnClean brush) is easy to use and designed to last up to 30 days. By using the MediLime ENClean brush or the MediLime disposable brushes, the ENFit connector remains clean, thereby feeding is more hygienic.

Article	Description	Quantity
 0901-03	MediLime EnClean brush	3
 ECB-NS	MediLime disposable brush	50

“  
Quality of life for all people.”



## Product range

MediLime offers a complete range of Tummy Buttons, Extension Tubes, Feeding Syringes and connectors. MediLime is continuously working on optimizing its range and strives to bring the best product for everyone.


Follow MediLime and keep up with the latest developments. [www.medilime.com](http://www.medilime.com)

# Medilime<sup>®</sup>

 Stavangerweg 21-14

 9723 JC Groningen

 The Netherlands

 +31 (0)85 40 19 648

 [info@medilime.com](mailto:info@medilime.com)

 [www.medilime.com](http://www.medilime.com)